

TASKI Sani 4 in 1 conc

Washroom cleaner, descaler, disinfectant and deodoriser

Description

Washroom cleaner, descaler, disinfectant and deodoriser for all acid resistant hard surfaces in washroom areas.

Key properties

- · Acid based formulation
- Cleans, descales, disinfects and deodorises in one step
- Effective against a wide spectrum of micro-organisms including E coli and Salmonella
- Unique, patented Odour Neutralizing Technology (O.N.T.)
- Part of the unique EXACT dosing and packaging concept

Benefits

- Excellent cleaner, works in a large spectrum of water hardness
- Reduces the number of products in the washroom area
- Excellent hygiene, reducing the risk of cross contamination
- Leaves a pleasant, fresh fragrance

Use instructions

Spray cleaning:

- 1. Spray onto a damp cloth/ sponge to apply product on surface and allow to act briefly (please refer to the microbiological tab).
- 2. Use a sponge to remove stubborn soil.
- 3. Rinse or replace the cloth regularly.

Bucket cleaning:

- 1. Apply solution with cloth/sponge or mop and wipe.
- 2. Let the solution act briefly (please refer to the microbiological tab).
- 3. Use a sponge to remove stubborn soil.

Dosage:

TASKI Sani 4 in 1 conc needs to be diluted.

Dosage spray: 1 cap/500ml. Dosage bucket: 4 caps/5L.

Increase dosage in case of heavy soiling.

Proper dosage and temperature save costs and minimize environmental impacts.









° A.I.S.



TASKI Sani 4 in 1 conc

Technical data

Appearance: Pink liquid pH value (neat): 0.2 pH value (in use): 1.2

Relative density (20°C): 1.044 g/cm³

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diversey.com. Store in original closed containers away from extremes of temperature. Only for professional users / specialists.

Product compatibility

Do not use on acid-sensitive surfaces, e.g. marble. Sufficiently pre-wet grout with water before use. Test material compatibility in a small, inconspicuous place before use.

Microbiological data

Bacteridial acc. to EN 1276 (1% / 5 min / dirty condition).

Bacteridial acc. to EN13697 (2% / 5 min / dirty condition).

Fungicidal acc. to EN13697 (2.5% / 15 min / dirty).

Yeasticidal acc. to EN13697 (1.5% / 15 min / dirty).

Virucidal acc. to EN14476 (Rotavirus) (1% / 10 min / clean condition).

Virucidal acc. to EN14476 (Influenza A virus H3N8) (1% / 1 min / dirty condition).

Virucidal acc. to EN14476 (Swine Influenza H1N1) ((1% / 1 min / dirty condition).

Yeasticidal acc. to EN1650 (0,5% / 15 min / dirty condition (only candida).

Fungicidal acc. to EN1650 (2.5% / 15 min / dirty).

Available pack sizes

TASKI Sani 4 in 1 conc is available in in 6x1L dosing bottle and 1x5L cubitainer.